

10-Position Autosampler

Available for:

Vapor Pressure Tester
ERAVAP

FTIR Fuel Analyzer
ERASPEC

Multi-Fuel Analyzer
ERASPEC OIL

Oil-in Water Tester
ERACHECK ECO

Oil in Water Tester
ERACHECK PRO

eralytics' **10-position autosampler** is an optional accessory for nearly all eralytics analyzers. The 10-position autosampler is directly attached to the instrument allowing to conveniently connect up to 10 samples simultaneously. It is the perfect solution for higher sample throughput and even more efficiency. All instruments maintain their portability, even with the 10-position autosampler attached.



FEATURES AT A GLANCE

- 10-position autosampler for high sample throughput
- True instrument portability, even with the 10-position autosampler attached
- Simultaneous connection of up to 10 samples
- All parameters (such as sample name, method, etc.) can be defined for each position individually
- Easiest modifying of pending sample positions during an ongoing measurement

SPECIFICATIONS

| | |
|------------------------|--|
| Sampling | 10 individual ports |
| Filling | Automatic via tubes or from glass syringes |
| Dimensions (W x H x D) | 24.2 x 15 x 7 cm (9.5 x 5.9 x 2.8 in) |
| Weight | 1.9 kg (4.1 lb) |