

# erasoft RCS Remote Control PC Software

The easy-to-use **erasoft** RCS software is the perfect tool to connect a PC running Microsoft Windows® with one or more **eralytics** analyzers. As an alternative to LIMS, it allows for convenient data transfer and remote control of the analyzers. The **erasoft** RCS (Standard and Premium Version) can be installed on multiple computers



**erasoft** RCS<sup>1</sup> is available in a Standard and Premium Version. The **Premium Version** offers many advanced features and unique functions, the **Standard Version** with basic functions is free of charge.

FEATURES	Standard	Premium
Comfortable result transfer to a PC	☑	☑
Optional interface to a company LIMS <sup>2</sup>	☑	☑
Convenient handling of output configuration files (printer, LIMS)	☑	☑
Log-files detailing all warnings and errors on the instrument	☑	☑
Live status display of the instrument	☑	☑
SpecView®: Graphically view and printing of spectra (for Eraspec and Eraspec Oil)		☑
Automatic result download		☑
Auto printing & LIMS export <sup>2</sup>		☑
Remote control for the measurement and corresponding parameters		☑
Comfortable remote software-update of the instrument		☑

<sup>1</sup> **erasoft** RCS is compatible with all **eralytics** instrument generations, but not all listed features may be available for older instrument generations.

<sup>2</sup> The RCS can be used as an interface if the direct connectivity of an instrument to a LIMS is not possible.